

Solid Hydrogen Target Cooling and Warming Procedure

December 2018

1 Introduction

This document contains the instructions for cooling and warming processes when forming the hydrogen target.

2 Setting up the Labview Program

1. Open Labview using SHT program shortcut on Desktop.
2. In Labview, change the logging file path. Define today's date in format yymmdd in place of xxxxxx and hit enter.
3. Click the arrow in the upper left corner of the window to start.
4. Click 'Configure' to send command to instrument, then click 'Run' to start running scan cycle for data logging.
5. In the upper right corner of the Labview window, set the scan interval (default 10 s) and the file-save interval (default 1/30).

The following list contains the labels and their descriptions from the Labview program.

VAC-CC10 Wide range vacuum monitor.

VAC-TM1 Hydrogen line vacuum monitor.

VAC-TM2 High vacuum monitor.

VAC-PM Low vacuum monitor.

PtCo1 Target temperature.

PtCo2 Radiation shield temperature.

TGT-Heater Target heater DC voltage.

H2-Press Hydrogen Pressure.

VAC-CC10 and **VAC-TM2** read the same vacuum using different gauges.

3 Cooling Process

Valve	Initial Condition	Purpose
V ₀	closed	Controls flow of H to diffuser.
V₁	closed	Fine control of H from system to V₀.
V₂	closed	Bypasses V₁ to hydrogen supply line and V₀.
V ₃	closed	Controls H flow from system to V ₂ .
V ₄	open	Controls flow from small hydrogen cylinder to system.
V ₅	closed	Opens system to roughing pump.
V ₆	closed	Controls flow of H ₂ /D ₂ from V ₇ to system.
V ₇	closed	Controls flow from H ₂ /D ₂ bottle to V ₆ .

Valves removed 2019